GAS

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1.A green triple-locked strategy based on volatile-compound imaging,chemometrics, and markers to discriminate winter honey and sapium honey using headspace gas chromatography-ion mobility spectrometry[J]. Food Research International, https://doi.org /10.1016/j.foodres. 2019.01.004.

2.Development of a flavor fingerprint by HS-GC–IMS with PCA for volatilecompounds of Tricholoma matsutake Singer. Food Chemistry, 290 (2019) 32-39.

3.Characterization of the aroma release and perception of white bread duringoral processing by gas chromatography-ion mobility spectrometry and temporal dominance of sensations analysis. Food Research International, 123 (2019) 612-622.

4.Robustness study of calibration models for olive oil classification: Targeted and non-targeted fingerprint approaches based on GC-IMS. Food Chemistry, 288 (2019) 315-324.

5.Direct authentication of three Chinese materia medica species of the Lilii Bulbus family in terms of volatile components by headspace-gaschromatography-ion mobility spectrometry. Analytical Methods, 2019. DOI: 10.1039/c8ay02338g.

6.Content variations in compositions and volatile component in jujube fruits during the blacking process. Food Sci Nutr. 2019. 1–9.

7.Identification of terpenes and essential oils by means of static headspace gas chromatography-ion mobility spectrometry. Anal Bioanal Chem, (2017)409 6595–6603.

8.Detection of Adulteration in Canola Oil by Using GC-IMS and Chemometric Analysis. International Journal of Analytical Chemistry, 2018, Article ID 3160265.

9.Characteristic volatiles fingerprints and changes of volatile compounds infresh and dried Tricholoma matsutake Singer by HS-GC-IMS and HS-SPME GC–MS. Journal of Chromatography B, 1099 (2018) 46-55.

10.Ion mobility spectrometry coupled to gas chromatography: A rapid tool toassess eggs freshness. Food Chemistry, 271 (2019) 691-696.

11.Volatile-Compound Fingerprinting by Headspace-Gas Chromatography Ion-Mobility Spectrometry (HS-GC-IMS) as a Benchtop Alternative to 1H NMR Profiling for Assessment of the Authenticity of Honey. Analytical Chemistry, DOI: 10.1021/acs.analchem.7b03748. 2018.

12.Coupling Gas Chromatography with Ion Mobility Spectrometry. LC*GC Europe, 2016. 294-303.

13.Target identification of volatile metabolites to allow the differentiation of lactic acid bacteria by gas chromatography-ion mobility spectrometry. Food Chemistry, 220 (2017) 362-370.

14.Ion mobility spectrometry of volatile compounds from Iberian pig fat for fast feeding regime authentication. Talanta, 76 (2008) 591–596.

15.A review of recent, unconventional applications of ion mobility spectrometry (IMS). Analytica Chimica Acta, 703 (2011) 114–123.

16.Using ion mobility spectrometry for screening the autoxidation of peanuts. Food Control, 64 (2016) 17-21.

17.Ion mobility spectrometry fingerprints: A rapid detection technology for adulteration of sesame oil. Food Chemistry, 192 (2016) 60-66.

18.Molecular Analysis for Determining Differences and Similarities of Perfume Dupes. GLOBAL INGREDIENTS &FORMULATIONS GUIDE 2016.

19.Ion Mobility Spectrometry versus Classical Physico-chemical Analysis for Assessing the Shelf Life of Extra Virgin Olive Oil According to Container Type and Storage Conditions. J. Agric. Food Chem, 63 (2015) 2179−2188.

20.Gas chromatography and associated techniques in the analysis of laundry malodour. H&PC Today - Household and Personal Care Today, Vol. 9 nr. 2014.

21.Evaluation and application of static headspace–multicapillary column-gas chromatography–ion mobility spectrometry for complex sample analysis. Journal of Chromatography A, 1338 (2014) 136–148.

22.Direct determination of 2,4,6-tricholoroanisole in wines by single-drop ionic liquid microextraction coupled with multicapillary column separation and ion mobility spectrometry detection. Journal of Chromatography A, DOI:10.1016/j.chroma.2011.06.032. 2011.
23.Volatile Metabolites of Goat Cheeses Determined by Ion Mobility Spectrometry. Potential Applications in Quality Control. Food Anal Methods, DOI: 10.1007/s12161-014-0050-1. 2014.
24.Multi-capillary column-ion mobility spectrometry: a potential screening system to differentiate virgin olive oils. Anal Bioanal Chem, 402 (2012) 489–498.
25.Ion Mobility Spectrometry a Versatile Analytical Tool for Metabolomics Applications in Food Science. Current Metabolomics, 2 (2014) 264-271.
26.Direct coupling of a gas–liquid separator to an ion mobility spectrometer for the classification of different white wines using chemometrics tools. Talanta, 84 (2011) 471–479.
27.Determination of volatile compounds by GC–IMS to assign the quality of virgin olive oil. Food Chemistry, 187 (2015) 572-579.
28.Resolution-optimized headspace gas chromatography-ion mobility spectrometry (HS-GC-IMS) for non-targeted olive oil profiling. Anal Bioanal Chem, DOI: 10.1007/s00216-017-0338-2. 2017.
29.Direct classification of olive oils by using two types of ion mobility spectrometers. Analytica Chimica Acta,DOI:10.1016/j.aca.2011.03.007.
30.Target vs spectral fingerprint data analysis of Iberian ham samples for avoiding labelling fraud using headspace – gas chromatography–ion mobility spectrometry.Food Chemistry, 246 (2018) 65-73.
31.Direct differentiation of herbal medicine for volatile components by a multicapillary column with ion mobility spectrometry method. J. Sep. Sci, 38 (2015) 3205–3208.
32.Use of a non-destructive sampling method for characterization of Iberian cured ham breed and feeding regime using GC-IMS. Meat Science, https://doi.org/10.1016/j.meatsci.2019.02.018. 2019.
33.Biochemical changes induced by dominant bacteria in chill-stored silver carp (Hypophthalmichthys molitrix) and GC-IMS identification of volatile organic compounds. Food Microbiology, 84 (2019) 103248.
34.HS-GC-IMS and chemometric data treatment for food authenticity assessment: Olive oil mapping and classification through two different devices as an example. Food Control, DOI: https://doi.org/10.1016/j.foodcont.2018.11.001. 2018.
35.Quality assessment of olive oils based on temperature-ramped HS-GC-IMS andsensory evaluation: Comparison of different processing approaches by LDA, kNN, and SVM. Food Chemistry, DOI: https://doi.org/10.1016/j.foodchem.2018.11.095. 2018.
36.Data fusion of GC-IMS data and FT-MIR spectra for the authentication of olive oils and honeys—is it worth to go the extra mile? Analytical and Bioanalytical Chemistry, https://doi.org/10.1007/s00216-019-01978-w, 2019.
37. Microbial volatile organic compounds and their application in microorganism identification in foodstuff. Trends in Analytical Chemistry. DOI: http://dx.doi.org/doi: 10.1016/j.trac.2015.08.010 2015.
38.Modification of gelatin hydrolysates from grass carp (Ctenopharyngodonidellus) scales by Maillard reaction: Antioxidant activity and volatile compounds. Food Chemistry, 295 (2019) 569-578.
39.Ion Mobility Spectrometry in Food Analysis: Principles, Current Applications and Future Trends. Molecules , 2019, 24, 2706; DOI :10.3390/molecules24152706.
40.Gas chromatography-ion mobility spectrometric classification of vegetable oils based on digital image processing. Journal of Food Measurement and Characterization, https://doi.org/10.1007/s11694-019-00116-5.

41.Evaluation of fast volatile analysis for detection of Botrytis cinerea infections in strawberry. Food Microbiology, 32 (2012) 406-414.
42.Rapid detection of hydrogen sulfide produced by pathogenic bacteria in focused growth media using SHS-MCC-GC-IMS. Microchemical Journal, 140 (2018) 232-240.
43.Variations of volatile flavour compounds in Cordyceps militaris chicken soup after enzymolysis pretreatment by SPME combined with GC-MS, GC ×GC-TOF -MS and GC-IMS. International Journal of Food Science and Technology. DOI:10.1111/ijfs.14294.
44.Coupling laser desorption with gas chromatography and ion mobility spectrometry for improved olive oil characterisation. Food Chemistry. 255 (2018) 323-331.
45.Changes in the Volatile Components of Candied Kumquats in Different Processing Methodologies with Headspace–Gas Chromatography–Ion Mobility Spectrometry. Molecules 2019, 24, 3053;DOI:10.3390/molecules24173053.
46.Gas Chromatography-Ion Mobility Spectrometry Detection of Odor Fingerprint as Markers of Rapeseed Oil Refined Grade. Journal of Analytical Methods in Chemistry. https://doi.org/10.1155/2019/3163204.

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1. Breath Analysis Using eNose and Ion Mobility Technology to Diagnose Inflammatory Bowel Disease—A Pilot Study. Biosensors 2019, 9, x; DOI: FOR PEER REVIEW.
2. Monitoring of selected skin- and breath-borne volatile organic compounds emitted from the human body using gas chromatography ion mobility spectrometry (GC-IMS). Journal of Chromatography B. 1076 (2018) 29-34.
3.Ion mobility spectrometry for microbial volatile organic compounds: a new identification tool for human pathogenic bacteria. Appl Microbiol Biotechnol.DOI: 10.1007/s00253-012-3924-4. 2012.
4.Detection and validation of volatile metabolic patterns over different strains of two human pathogenic bacteria during their growth in a complex medium using multi-capillary column-ion mobility spectrometry (MCC-IMS). Appl Microbiol Biotechnol. DOI: 10.1007/s00253-013-4762-8. 2013.
5.Analysis of Listeria using exogenous volatile organic compound metabolites and their detection by static headspace–multi-capillary column–gas chromatography–ion mobility spectrometry (SHS–MCC–GC–IMS). Anal Bioanal Chem. DOI: 10.1007/s00216-017-0375-x. 2017.
6.Comparison of metabolites in exhaled breath and bronchoalveolar lavage fluid samples in a mouse model of asthma. J Appl Physiol. 111(2011) 1088–1095.
7.Detection of human metabolites using multi-capillary columns coupled to ion mobility spectrometry. Journal of Chromatography A. 1084 (2005) 145-151.
8.Detection of Metabolites of Trapped Humans Using Ion Mobility Spectrometry Coupled with Gas Chromatography. Anal. Chem. dx.doi.org/10.1021/ac302752f . 2012.
9.Determination of Volatile Organic Compounds in Exhaled Breath by Ion Mobility Spectrometry. Chem. Anal. (Warsaw), 53 (2008) 953.

10.Ion Mobility Spectrometry for Detection of Skin Volatiles. Journal of Chromatography B. DOI:10.1016/j.jchromb.2012.10.028. 2012.
11.Ion mobility spectrometry for pharmacokinetic studies—exemplary application. J. Breath Res. 7 (2013) 046008 (9pp) DOI: 10.1088/1752-7155/7/4/046008.
12.Potential and suitability of Ion Mobility Spectrometry (IMS) for breath analysis. Technology and Medical Sciences – Natal Jorge et al. (eds).2011.
13.Volatile organic compounds – a new class of uremic toxins.
14.Non-invasive monitoring of bacterial growth and auto-induced protein production in a bioreactor with a closed-loop GC-IMS. Int. J. Ion Mobil. Spec. 18 (2015) 9–15. DOI: 10.1007/s12127-014-0163-7.
15.Ananysis of human breath using IMS. IJIMS 8 (2005) 1, 5-8.
16.Determination of microbial volatile organic compounds (MVOC)using IMS with different ionisation sources. IJIMS. 5 (2002) 3, 138-142, p. 139.
17.Peak assignment in multi-capillary column–ion mobility spectrometry using comparative studies with gas chromatography–mass spectrometry for VOC analysis. Anal Bioanal Chem. 396 (2010) 471–482.
18.Detection of characteristic metabolites of Aspergillus fumigatus and Candida species using ion mobility spectrometry – metabolic profiling by volatile organic compounds. Mycoses 54, e828–e837.2011.
19.Direct authentication of three Chinese materia medica species of the Lilii Bulbus family in terms of volatile components by headspace-gas chromatography-ion mobility spectrometry. Analytical Method.,2019; DOI: 10.1039/c8ay02338g.
20.Direct differentiation of herbal medicine for volatile components by a multicapillary column with ion mobility spectrometry method. J. Sep. Sci. 38(2015) 3205–3208.
21.Measurement of exhaled volatile organic compounds from patients with chronic obstructive pulmonary disease (COPD) using closed gas loop GC-IMS and GC-APCI-MS. J. Breath Res. 2016,10, 026004.
22.Discrimination of bacteria by rapid sensing their metabolic volatiles using an aspiration-type ion mobility spectrometer (a-IMS) and gas chromatography-mass spectrometry GC-MS. Analytica Chimica Acta. 2017, 1-9.
23.A volatile biomarker in breath predicts lung cancer and pulmonary nodules. Journal of Breath Research. 2019. https://doi.org/10.1088/1752-7163/ab21aa.
24.A novel chromatographic approach to distinguish Gram-positive from Gram-negative bacteria using exogenous volatile organic compound metabolites. Journal of Chromatography A. http://dx.doi.org/doi:10.1016/j.chroma.2017.04.015. 2017.
25.Extraction of toxic compounds from saliva by magnetic-stirring-assisted micro-solid-phase extraction step followed by headspace-gas chromatography-ion mobility spectrometry. Anal Bioanal Chem. (2016) 408:6813–6822.
26.Trace detection of endogenous human volatile organic compounds for search, rescue and emergency applications. Trends in Analytical Chemistry. 66 (2015) 158-175.
27.Bedside Measurement of VolatileOrganic Compounds in theAtmosphere of Neonatal IncubatorsUsing Ion Mobility Spectrometry. Frontiers in Pediatrics. DOI: 10.3389/fped.2019.002.

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Reference: Rapid discrimination of Citrus reticulata ‘Chachi’ by headspace-gas chromatography-ion mobility spectrometry fingerprints combined with principal component analysis. **Food Research International**. 131 (2020) 108985.

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
Fig. 5. Topographic plot of GC–IMS spectra. (A) Pomace olive oil, (B) olive oil and (C) extra virgin olive oil.


Fig. 3. PCA results for mean-centered data for the three types of oil as obtained with UV–IMS. Samples: Pomace olive oil (red triangles); olive oil (green asterisks); extra virgin olive oil (blue squares). (For interpretation of the references to color in this figure legend, the reader is referred to the web version of the article.)

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Reference: Direct classifification of olive oils by using two types of ion mobility spectrometers. Anal. Chim.**Acta (2011).DOI: 10.1016/j.aca.2011.03.007.**

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Referenc：Volatile-Compound Fingerprinting by Headspace-Gas-Chromatography Ion-Mobility Spectrometry (HS-GC-IMS) as aBenchtop Alternative to 1H NMR Profiling for Assessment of the **Authenticity of honey. Analytical Chemistry DOI: 10.1021/acs.analchem.7b03748（2017）**

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Reference：Ion mobility spectrometry coupled to gas chromatography: A rapid tool toassess eggs freshness. **Food Chemistry, 271 (2019) 691-696.**

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Reference: Characteristic volatiles fingerprints and changes of volatile compounds in fresh and dried Tricholoma matsutake Singer by HS-GC-IMS and HS-SPME-GC–MS. **Journal of Chromatography B, 1099. (2018) 46-55.**

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Reference：Content variations in compositions and volatile component in jujube fruits during the blacking process. **Food Sci Nutr. 2019;1-9.**

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8. Discrimination of Chinese Yellow Wine from Different Origins Based on Flavor Fingerprint. **Acta Chromatographica**. DOI: 10.1556/1326.2019.00613.

9. Rapid discrimination of Citrus reticulata 'Chachi' by headspace-gas chromatography-ion mobility spectrometry fingerprints combined with principal component analysis. **Food Research International**. 131 (2020) 108985. https://doi.org/10.1016/j.foodres.2020.108985.

10. Stable isotope signatures versus gas chromatography ion mobility spectrometry to determine the geographical origin of Fujian Oolong tea (Camellia sinensis) samples. **European Food Research and Technology** (2020) 246:955–964. https://doi.org/10.1007/s00217-020-03469-0

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5. Direct authentication of three Chinese materia medica species of the Lilii Bulbus family in terms of volatile components by headspace-gaschromatography-ion mobility spectrometry. **Analytical Methods**. DOI: 10.1039/c8ay02338g.

6. Volatile-Compound Fingerprinting by Headspace Gas Chromatography Ion-Mobility Spectrometry (HS-GC-IMS) as a Benchtop Alternative to 1H NMR Profiling for Assessment of the Authenticity of Honey. **Anal. Chem**. DOI: 10.1021/acs.analchem.7b03748.

7. Untargeted and Targeted Discrimination of Honey Collected by Apis cerana and Apis mellifera Based on Volatiles Using HS-GC-IMS and HS-SPME-GC-MS. Journal of Agricultural and Food Chemistry. **J. Agric. Food Chem**. 2019, 67, 12144—12152.

8. Identification of terpenes and essential oils by means of static headspace gas chromatography-ion mobility spectrometry. **Anal Bioanal Chem** (2017) 409:6595–6603.

9. Volatile Metabolites of Goat Cheeses Determined by Ion Mobility Spectrometry Potential Applications in Quality Control. **Food Anal. Methods**. DOI 10.1007/s12161-014-0050-1.

10. Characterization of Arabica and Robusta Coffees by Ion Mobility Sum Spectrum. **Sensors** 2020, 20, 3123; doi:10.3390/s20113123.

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5. Direct differentiation of herbal medicine for volatile components by a multicapillary column with ion mobility spectrometry method. **J. Sep. Sci**. 2015, 38, 3205–3208.

6. Target vs spectral fingerprint data analysis of Iberian ham samples for avoiding labelling fraud using headspace – gas chromatography–ion mobility spectrometry. **Food Chemistry**. 246 (2018) 65-73.

7. A green triple-locked strategy based on volatile-compound imaging,chemometrics, and markers to discriminate winter honey and sapium honey using headspace gas chromatography-ion mobility spectrometry. **Food Research International**. https://doi.org/10.1016/j.foodres.2019.01.004.

8. Untargeted headspace gas chromatography – Ion mobility spectrometry analysis for detection of adulterated honey. **Talanta**. 205 (2019) 120123.

9. Detection of Peanut Oil Adulteration Mixed with Rapeseed Oil Using Gas Chromatography and Gas Chromatography–Ion Mobility Spectrometry. **Food Analytical Methods**.
https://doi.org/10.1007/s12161-019-01571-y.

10. Molecular Analysis for Determining Differences and Similarities of Perfume Dupes. **Global Ingedients &Formulations Gulde**. 2016.

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4. A robustness study of calibration models for olive oil classification: Targeted and non-targeted fingerprint approaches based on GC-IMS. **Food Chemistry**. 288 (2019) 315-324.

5. Resolution-optimized headspace gas chromatography-ion mobility spectrometry (HS-GC-IMS) for non-targeted olive oil profiling. **Anal Bioanal Chem**. DOI 10.1007/s00216-017-0338-2.

6. Direct classification of olive oils by using two types of ion mobility spectrometers. **Anal. Chim. Acta** (2011), doi:10.1016/j.aca.2011.03.007.

7. Gas Chromatography-Ion Mobility Spectrometry Detection of Odor Fingerprint as Markers of Rapeseed Oil Refined Grade. **Journal of Analytical Methods in Chemistry**. https://doi.org/10.1155/2019/3163204.

8. Usefulness of GC-IMS for rapid quantitative analysis without sample treatment: Focus on ethanol, one of the potential classification markers of olive oils. **LWT - Food Science and Technology**. 120 (2020) 108897.

9. Use of headspace-gas chromatography-ion mobility spectrometry to detect volatile fingerprints of palm fibre oil and sludge palm oil in samples of crude palm oil. **BMC Res Notes** (2019) 12:229 https://doi.org/10.1186/s13104-019-4263-7.

10. Characteristic volatiles fingerprints and profiles determination in different grades of coconut oil by HS-GC-IMS and HS-SPME-GC-MS. **International Journal of Food Science and Technology**. 2020. doi:10.1111/ijfs.14664.

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5. Comparison of fermentation behaviors and properties of raspberry wines by spontaneous and controlled alcoholic fermentations. **Food Research International**. 128 (2020) 108801.

6. Changes in the Volatile Components of Candied Kumquats in Different Processing Methodologies with Headspace–Gas Chromatography–Ion Mobility Spectrometry. **Molecules** 2019, 24, 3053; doi:10.3390/molecules24173053.

7. Evaluation by electronic tongue and headspace-GC-IMS analyses of the flavor compounds in dry-cured pork with different salt content. **Food Research International**. https://doi.org/10.1016/j.foodres.2020.109456.

8. Biochemical changes induced by dominant bacteria in chill-stored silver carp (Hypophthalmichthys molitrix) and GC-IMS identification of volatile organic compounds. **Food Microbiology**. 84 (2019) 103248.

9. Mechanism of aroma compounds changes from sea cucumber peptide powders (SCPPs) under different storage conditions. **Food Research International**. 128 (2020) 108757.

10. Changes in volatile flavor compounds of peppers during hot air drying process based on headspace-gas chromatography-ion mobility spectrometry (HS-GC-IMS). **J Sci Food Agric**. 2020; 100: 3087–3098.

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4. Correlation analysis of the age of brandy and volatiles in brandy by gas chromatography-mass spectrometry and gas chromatography-ion mobility spectrometry. **Microchemical Journal**. 157 (2020) 104948.

5. Characterization of Volatile Component Changes in Jujube Fruits during Cold Storage by Using Headspace-Gas Chromatography-Ion Mobility Spectrometry. **Molecules** 2019, 24, 3904; doi:10.3390/molecules24213904.

6. Volatile Metabolites of Goat Cheeses Determined by Ion Mobility Spectrometry. Potential Applications in Quality Control. **Food Anal Methods**. DOI 10.1007/s12161-014-0050-1.7

7. Characterization of the aroma release and perception of white bread during oral processing by gas chromatography-ion mobility spectrometry and temporal dominance of sensations analysis. **Food Research International**. 123 (2019) 612-622.

8. Ion mobility spectrometry coupled to gas chromatography: A rapid tool to assess eggs freshness. **Food Chemistry**. 271 (2019) 691-696.
9. Analysis of the volatile compounds associated with pickling of ginger using headspace gas chromatography – ion mobility spectrometry. **Flavour Fragr J**. DOI: 10.1002/ffj.3530.
10. Ion Mobility Spectrometry versus Classical Physico-chemical Analysis for Assessing the Shelf Life of Extra Virgin Olive Oil According to Container Type and Storage Conditions. **J. Agric. Food Chem**. 2015, 63, 2179—2188.

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1. 臭鳜鱼低温发酵过程中品质及挥发性物质的变化. **中国食品学报**. 2020.
2. 气相离子迁移谱联用技术评定大米霉变程度的研究应用. **中国粮油学报**. 2019.
3. Characterization of aroma profile and characteristic aromas during lychee wine fermentation. **J Food Process Preserv**. DOI: 10.1111/jfpp.14003.
4. A Study on Volatiles Metabolites Screening by HS-SPME-GC-MS and HS-GC-IMS for Discrimination and Characterization of White and Yellowed Rice. **Cereal Chemistry**. 2020; 00:1–9.. DOI: 10.1002/cche.10264.
5. Evaluation of fast volatile analysis for detection of Botrytis cinerea infections in strawberry. **Food Microbiology** 32 (2012) 406e414.
6. Effects of different probiotic combinations on the components and bioactivity of Spirulina. **J Basic Microbiol**. 2020;1–15.
7. Changes in volatile compounds of fermented minced pepper during natural and inoculated fermentation process based on headspace–gas chromatography–ion mobility spectrometry. **Food Sci Nutr**. 2020; 00:1–18.
8. Target identification of volatile metabolites to allow the differentiation of lactic acid bacteria by gas chromatography-ion mobility spectrometry. **Food Chemistry**. 220 (2017) 362–370.
9. Rapid detection of hydrogen sulfide produced by pathogenic bacteria in focused growth media using SHS-MCC-GC-IMS. **Microchemical Journal**. 140 (2018) 232–240.
10. Profiling and characterization of odorous volatile compounds from the industrial fermentation of erythromycin. **Environmental Pollution**. 255 (2019) 113130.

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